

Operating instructions
Istruzioni per l'uso
Bedienungsanleitung
Mode d'emploi
Modo de empleo

REF 234.52

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MULTiflex® is a registered trademark of Kaltenbach & Voigt GmbH, Germany.

TeKne Dental s.r.l.

Via del Pescinale, 77 - 50041 Calenzano (FI) - Italy
☎ +39 055 8825741 📠 +39 055 8825764
info@teknedental.com www.teknedental.com

DESCRIPTION

Pneumatic handpiece, compliant to ISO 15606 Standard.
The handpiece can be used for prophylaxis procedures or bone surgery employing the Komet Sonosurgery Inserts.
Weight without insert: 88 g
Frequency of vibrations: 5000 ÷ 6500 Hz.
Maximum noise level: 75 dBA (@3.6 bar)

This medical device meets the requirements of the European Directives 93/42/EEC (Class IIa) and 2007/47/EC.

CONTENTS

SONOSURGERY® handpiece, a insert-fitting dynamometric wrench.

CONNECTION

The SONOSURGERY® air handpiece can be connected to any unit supplied with a standard turbine hose.
The SONOSURGERY® handpiece with quick connection (Fig. 2) must be used together with the SONOSURGERY® STERIL (REF 316.40.08) or MULTiflex® quick couplings.

AIR PRESSURE

Recommended air pressure is 2.8 ÷ 3.6 bar (40.5 ÷ 52.1 psi).
Maximum air consumption is 22 NI/min (@3.6 bar).
Check air pressure by means of a gauge (REF 129.40) connected between the handpiece and the supply hose.
Supplied air must be perfectly dry and well-filtered, kept at constant pressure.
Impurity or residues in the compressed air may cause the handpiece to malfunction.

COOLING

Operating water pressure for cooling spray must be 0.5 ÷ 2.5 bar. Water flow rate is higher than 50 ml/min (@2.0 bar). Water must be well-filtered.

FITTING OF INSERTS

Select the right type of insert according to the operation to be performed and screw it manually onto the handpiece. The insert must then be moderately tightened using the supplied key (Fig. 1).
Inserts must always be removed after use.
Do not fit inserts while the handpiece is running.

POWER ADJUSTMENT

Power may be increased or decreased by rotating the ring nut located at the bottom of the handpiece (Fig. 3).

MAINTENANCE OF INSERTS

The shape and weight of each insert are the result of a precise design optimisation so as to obtain maximum performance.
Therefore any alteration, breaking, twisting or wearing of the insert can cause abnormal operation: in such case it is strongly recommended that the insert be replaced.

LUBRICATION

In order for the instrument you have purchased to be helpful in your profession for a very long time, carefully repeat the following operation **weekly**.
Detach the handpiece from the supply hose and lubricate using the special can (REF 120.00). The SO2000 service oil must be inserted into instrument air input. Press once only.
In case of handpiece with quick connection, insert the special nozzle reducer (REF 120.02) and lubricate as shown on Fig. 4.
Immediately after, screw the handpiece back to the hose (or reconnect it to the quick coupling) and operate the instrument for a few seconds, until the characteristic hissing sound can be heard.

CLEANING

Manual external: Carefully clean the surfaces using a clean cloth soaked in a suitable substance. Use also a medium toothbrush and brush off under flowing hot tap water.
Manual internal: Using a suitable spray product, follow corresponding instruc-tions and spray the substance to the internal of the medical device. Immediately after internal cleaning, carry out the dis-infection process or the drying process.
Automated: Use a thermodisinfecter complying with ISO 15883 standard which is operated with alkaline cleaning agents at a pH of maximum 10.

DISINFECTION

Manual external: Carefully clean the surfaces using a clean cloth soave in a suitable substance (i.e. O-phenylphenol or alcohols based on ethanol). Products not recommended: disinfectants con-taining benzalkonium chloride, acetone or glutaraldehyde.
Manual internal: Using a suitable spray product, follow corresponding instruc-tions and spray the substance to the internal of the medical device. Immediately after internal disinfection, carry out the OIL lubrication process and/or the sterilization process.
Automated: Use a thermodisinfecter complying with ISO 15883 standard which is operated with alkaline disinfectant agents at a pH of maximum 10.

DRYING

Manual: In general use dry and clean air complying to ISO 7494-2 standard.
Automated: The drying phase is normally part of the cleaning program of the thermodisinfecter.

STERILIZATION

 The handpiece, the wrench and inserts can be sterilized in autoclave at 134 °C/2 bar (273 °F/29 psi).
The handpiece must always be sterilized without the insert on.
Before sterilization, remove residues and carefully clean the instrument.

 During sterilization, metal parts of different kinds should not come into contact with each other. To avoid this possibility, wrap the objects in a cloth or place them in a sterilisable bag.

The SONOSURGERY® STERIL quick couplings can be sterilized.

GARANTEE

The device carry a 12-month guarantee against all defects of construction. Guar-antee validity is subject to specific conditions.

SERVICE

In case of malfunction and for any overhaul and repair work, please contact your usual supplier or apply directly to TKD.

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